
DE Gebrauchsanweisung **EN** Instructions for Use **FR** Mode d'emploi
NL Gebruiksaanwijzing **IT** Istruzioni d'uso **TR** Kullanma Kilavuzu
ES Manual de instrucciones


WM 68290a 05/2019 DE, EN, FR, NL, IT, TR, ES (MX)

CARA Full Face

Full Face Mask








**Löwenstein Medical
Technology GmbH + Co. KG**
Kronsalsweg 40
22525 Hamburg, Germany
T: +49 40 54702-0
F: +49 40 54702-461
www.loewensteinmedical.com

CE 0197

WM 68290a

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Product	Warranty periods
Masks including accessories	6 months

12 Declaration of conformity

The manufacturer Löwenstein Medical Technology GmbH + Co. KG (KronsaaIsweg 40, 22525 Hamburg, Germany) hereby declares that the product complies with the relevant regulations of Directive 93/42/EEC governing medical devices. The un-abridged text of the Declaration of Conformity can be found on the manufacturer's website.

1 Operation

The following figures show you how to adjust, apply, remove, dismantle and assemble the mask:

- 1** Putting on the mask
- 2** Adjusting the mask
- 3** Removing the mask
- 4** Dismantling the mask
- 5** Assembling the mask

2 Introduction

2.1 Intended use

The CARA Full Face mask is used for treating sleep apnea and for non-invasive and non-life-sustaining ventilation of patients with respiratory insufficiency. It serves as a connecting element between the patient and the therapy device.

2.2 Contraindications

The mask may not be used in the following situations: necessity for immediate intubation; loss of consciousness, acute vomiting.

The mask may be used in the following situations only with particular caution: pressure points and acute injuries to the facial skin; skin allergies involving the face; deformities of the face or nasopharynx; acute pain affecting the face; cough reflex restricted or absent; claustrophobia; acute nausea.

If you are not certain whether one of these situations applies to you, please consult your attending physician or medical advisor. Please also observe the contraindications in the instructions for use of your therapy device.

2.3 Side effects

The following side effects may occur with use of the mask: nasal congestion, dry nose, dry mouth in the morning, feeling of pressure in the sinuses, irritated conjunctiva, skin rashes, pressure marks on the face, irritating noises when breathing. Should such side effects occur, please contact your attending physician or medical advisor.

3 Safety

Risk of injury from re-inhaling CO₂

If the mask is used incorrectly, CO₂ may be re-inhaled.

- ⇒ Do not close off the exhalation system of the mask.

- ⇒ Only put on the mask for an extended period if the therapy device is running.
- ⇒ Only use the mask within the quoted therapy pressure range.
- ⇒ Patients who are unable to take the mask off themselves must be monitored by qualified nursing staff.

Risk of injury if the mask slips!

If the mask slips or falls off, the therapy is ineffective.

- ⇒ Monitor patients with restricted spontaneous respiration.
- ⇒ Activate low pressure/leakage alarms on the therapy device.

Risk of injury from anesthetic gases!

Anesthetic gas may escape through the exhalation valve and put third parties at risk.

- ⇒ Never use the mask during anesthesia.

Risk of injury from lack of cleaning!

The mask may show contamination, which can possibly put the patient at risk.

- ⇒ Clean the mask before using for the first time (see section entitled "Hygiene treatment").
- ⇒ Clean the mask regularly.

4 Product description

4.1 Overview

A diagram of the individual parts can be found on the title page.

1. Headgear
2. Link
3. Elbow
4. Anti-asphyxia valve
5. Rotating sleeve
6. Mask body
7. Headgear clip
8. Mask cushion

4.2 Compatible devices

With many combinations of device, actual pressure in the mask will not correspond to prescribed therapy pressure, even if the therapy device is displaying the correct pressure. Have the combination of devices adjusted by a physician or specialist dealer so that actual pressure in the mask corresponds to therapy pressure.

4.3 Exhalation system

The mask has an integrated exhalation system. The link and mask body are shaped so that there is a gap

between these parts. The exhaled air can escape through this gap.

4.4 Anti-asphyxia valve

▲ WARNING

Risk of asphyxia if anti-asphyxia valve not working properly!

Residues may cause the valve to stick and lead to CO₂ being re-inhaled.
⇒ Check before every use that the openings of the anti-asphyxia valve are clear.

if the therapy device fails, the anti-asphyxia valve opens so that the patient can breathe ambient air.

5 Hygiene treatment

▲ WARNING

Risk of injury from inadequate cleaning!

Residues may congest the mask, impair the integrated exhalation system and jeopardize therapy success.

⇒ For patients with a compromised immune system or particular background of illness, disinfect mask components daily following consultation with the physician.

5.1 Clean mask

1. Dismantle mask (see Figure 4).
2. Clean mask in accordance with the table below:

Action	Daily	Weekly
Wash mask components with warm water and a mild detergent.	X	
When washing mask parts, clean thoroughly with a cloth or a soft brush. Or: Put the mask parts in the dishwasher. Wash headgear by hand.		X
3. Rinse all parts with clear water.		X
4. Allow all parts to air-dry.		
5. Perform a visual inspection.		
6. If necessary: replace damaged parts.		
7. Re-assemble mask (see Figure 5).		

i Discolorations of mask components do not impair the functionality of the mask.

5.2 Change of patient

In the clinical sphere only: You must subject the mask to a hygiene treatment in the event of a change of patient. Information on the hygiene treatment with a change of patient can be found in a brochure available on the manufacturer's website. We will also send you this brochure on request.

6 Disposal

You can dispose of all parts in domestic waste.

7 Troubleshooting

Fault	Cause	Remedy
Pain from pressure on the face.	Mask too tight.	Loosen headgear slightly.
Draft in the eye.	Mask too loose. Mask does not fit.	Tighten headgear slightly. Contact your specialist dealer.
	Mask not correctly adjusted.	Re-adjust mask (see Figure 2).
	Mask cushion damaged.	Replace mask cushion.
Therapy pressure is not reached.	Patient circuit leaking. Anti-asphyxia valve defective (only vented).	Check connectors, check tubes properly located. Replace the mask.

8 Technical specifications


Product class to directive 93/42/EEC	IIa
Dimensions (W x H x D)	93 mm x 145 mm x 89 mm 95 mm x 160 mm x 90 mm 95 mm x 174 mm x 91 mm
Size S	
Size M	
Size L	
Weight	93 g 97 g 102 g
Dead space	180 ml 219 ml 244 ml
Size S	
Size M	
Size L	
Therapy pressure range	4 hPa - 25 hPa
Tube connection: tapered connection to EN ISO 5356-1 vented	Ø 22 mm (male)
Temperature range:	+5 °C to +40 °C -20 °C to +70 °C
Operation	
Transport and storage	
Flow resistance at 50 l/min	0.15 hPa
at 100 l/min	0.5 hPa
Flow resistance, anti-asphyxia valve	
Inspiration at 50 l/min:	0.6 hPa
Exhalation at 50 l/min:	0.8 hPa
Switching pressure	
Anti-asphyxia valve	
• Open:	0.5 hPa
• Close:	2.2 hPa
Quoted two-figure noise emission value to ISO 4871:	
- Sound pressure level	19 dB(A)
- Sound power level	27 dB(A)
- Uncertainty factor	3 dB(A)
Service life	up to 12 months ¹
Standards applied	EN ISO 17510-2:2009

¹ The mask materials deteriorate if exposed to e.g. aggressive detergents. In individual cases it may be necessary to replace mask parts sooner.

The right to make design modifications is reserved.

9 Pressure/flow curve

The characteristic pressure/flow curve shows the leakage flow as a function of therapy pressure.

10 Materials

Mask part	Material
Mask cushion	SI (silicone)
Link	PA (polyamide)
Headgear clip	PA (polyamide)
Rotating sleeve, mask body, elbow	PA (polyamide)
Headgear	Elastane, polyester, PU (polyurethane), CO (cotton)
Emergency exhalation system: Elbow, anti-asphyxia valve, valve safety device	PA (polyamide), SI (silicone), PP (polypropylene)

No parts of the mask contain latex, PVC (polyvinyl chloride) or DEHP (diethylhexyl phthalate).

11 Warranty

Löwenstein Medical gives the customer a limited manufacturer warranty on new original Löwenstein Medical products and any replacement part fitted by Löwenstein Medical in accordance with the warranty conditions applicable to the product in question and in accordance with the warranty periods from date of purchase as listed below. The warranty conditions are available on the website of the manufacturer. We will also send you the warranty conditions on request.

Please bear in mind that any claim to warranty and liability shall be void if neither the accessories recommended in the instructions for use nor genuine replacement parts are used.

In the event of a claim under warranty, contact your specialist dealer.

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